Safety Communication

Recall of Laser photocoagulator Easyret manufactured
by Quantel Medical.

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Device Name:

Laser photocoagulator Easyret with software version
1.03.

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Problem:

A software bug lead to delivery of an erroneous level of
energy corresponding to a 100% duty cycle when
treatment session be created. This issue could cause
photocoagulation of the retina in patients.

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Recommendation:

- Discontinue use of the MicroPulse mode on the
device until the software update in version 1.04
has been completed.
- Contact the manufacturer's authorized
representative in KSA on contact info below to
update the software.
- Share this information with Healthcare
Professionals in your facility.

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The affected devices are:

Laser photocoagulator Easyret
Software version 1.03
Product Code: QM-PABX0017
Serial Number: 155

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You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative for corrective action.

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Healthcare Professionals should **report** any adverse events suspected to be associated with affected devices above (or other Medical Devices) to:

National Center for Medical Devices Reporting (NCMDR Website).

Medical Devices Sector
Saudi Food and Drug Authority
Postal Address: Saudi Arabia - Saudi Food and Drug Authority (3292)
North Ring Road - Al Nafal Unit (1)
Riyadh 13312 - 6288
Tel: +966 (11) 2038222 Ext: 2406
Fax: +966 (11) 2757245

For latest published Recalls/Alerts, please visit NCMDR website:

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